

Forkhead box protein M1

Catalog No: #AP83473



Package Size: #AP83473-1 50ug #AP83473-2 100ug #AP83473-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Forkhead box protein M1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Forkhead-related protein FKHL16,Hepatocyte nuclear factor 3 forkhead homolog 11,HFH-11,M-phase phosphoprotein 2,MPM-2 reactive phosphoprotein 2,Transcription factor Trident,Winged-helix factor from INS-1 cells,FKHL16,HFH11,MPP2,WIN
Accession No.	Uniprot ID:Q08050
Uniprot	Q08050
GeneID	2305;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FOXM1

Note: This product is for in vitro research use only