

Heat shock-related 70 kDa protein 2

Catalog No: #AP83477



Package Size: #AP83477-1 50ug #AP83477-2 100ug #AP83477-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Heat shock-related 70 kDa protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Heat shock protein 70.2,Hcp70.2,Hsp70-2
Accession No.	Uniprot ID:P17156
Uniprot	P17156
GeneID	15512;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Hspa2

Note: This product is for in vitro research use only