

Golgi phosphoprotein 3

Catalog No: #AP83481



Package Size: #AP83481-1 50ug #AP83481-2 100ug #AP83481-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Golgi phosphoprotein 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Coat protein GPP34,Mitochondrial DNA absence factor,MIDAS,GPP34
Accession No.	Uniprot ID:Q9H4A6
Uniprot	Q9H4A6
GeneID	64083;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:GOLPH3

Note: This product is for in vitro research use only