

Galanin peptides

Catalog No: #AP83486



Package Size: #AP83486-1 50ug #AP83486-2 100ug #AP83486-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Galanin peptides
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	GAL1,GALN,GLNN
Accession No.	Uniprot ID:P22466
Uniprot	P22466
GeneID	51083;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:GAL

Note: This product is for in vitro research use only