

Histone H3.1t

Catalog No: #AP83489



Package Size: #AP83489-1 50ug #AP83489-2 100ug #AP83489-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Histone H3.1t
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	H3/t,H3/g,H3FT
Accession No.	Uniprot ID:Q16695
Uniprot	Q16695
GeneID	8290;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:HIST3H3

Note: This product is for in vitro research use only