

Proto-oncogene c-Fos

Catalog No: #AP83492



Package Size: #AP83492-1 50ug #AP83492-2 100ug #AP83492-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Proto-oncogene c-Fos
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	Cellular oncogene fos
Accession No.	Uniprot ID:P11939
Uniprot	P11939
GeneID	396512;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: FOS

Note: This product is for in vitro research use only