

Glucosylceramidase

Catalog No: #AP83493



Package Size: #AP83493-1 50ug #AP83493-2 100ug #AP83493-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Glucosylceramidase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Acid beta-glucosidase,Beta-glucocerebrosidase,D-glucosyl-N-acylsphingosine glucohydrolase
Accession No.	Uniprot ID:Q2KHZ8
Uniprot	Q2KHZ8
GeneID	537087;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:GBA

Note: This product is for in vitro research use only