

# Glutathione S-transferase A4

Catalog No: #AP83500



Package Size: #AP83500-1 50ug #AP83500-2 100ug #AP83500-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Glutathione S-transferase A4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	GST A4-4,GSTA4-4,GST class-alpha member 4,Glutathione S-transferase 5.7,GST 5.7,Gsta
Accession No.	Uniprot ID:P24472
Uniprot	P24472
GeneID	14860;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Gsta4

Note: This product is for in vitro research use only