

G-protein coupled receptor family C group 6 member A

Catalog No: #AP83502

Package Size: #AP83502-1 50ug #AP83502-2 100ug #AP83502-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	G-protein coupled receptor family C group 6 member A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Accession No.	Uniprot ID:E1BPQ3
Uniprot	E1BPQ3
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: GPRC6A

Note: This product is for in vitro research use only