

Isocitrate dehydrogenase [NADP] cytoplasmic

Catalog No: #AP83507



Package Size: #AP83507-1 50ug #AP83507-2 100ug #AP83507-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Isocitrate dehydrogenase [NADP] cytoplasmic
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	IDH, Cytosolic NADP-isocitrate dehydrogenase, IDP, NADP(+)-specific ICDH, Oxalosuccinate decarboxylase, PICD
Accession No.	Uniprot ID:O75874
Uniprot	O75874
GeneID	3417;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: IDH1

Note: This product is for in vitro research use only