

Furin

Catalog No: #AP83512



Package Size: #AP83512-1 50ug #AP83512-2 100ug #AP83512-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Furin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Dibasic-processing enzyme, Paired basic amino acid residue-cleaving enzyme, PACE, FUR, PACE, PCSK3
Accession No.	Uniprot ID:P09958
Uniprot	P09958
GeneID	5045;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: FURIN

Note: This product is for in vitro research use only