

Frizzled-8

Catalog No: #AP83513



Package Size: #AP83513-1 50ug #AP83513-2 100ug #AP83513-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Frizzled-8
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Fz-8
Accession No.	Uniprot ID:Q61091
Uniprot	Q61091
GeneID	14370;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Fzd8

Note: This product is for in vitro research use only