

# Pancreatic secretory granule membrane major glycoprotein GP2

Catalog No: #AP83516

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP83516-1 50ug #AP83516-2 100ug #AP83516-3 1mg

## Description

Product Name	Pancreatic secretory granule membrane major glycoprotein GP2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Glycoprotein 80, gp80, Pancreatic zymogen granule membrane protein GP-2
Accession No.	Uniprot ID:P19218
Uniprot	P19218
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Gp2

Note: This product is for in vitro research use only