

Granzyme B

Catalog No: #AP83517



Package Size: #AP83517-1 50ug #AP83517-2 100ug #AP83517-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Granzyme B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	C11,CTLA-1,Cathepsin G-like 1,CTSLG1,Cytotoxic T-lymphocyte proteinase 2,Lymphocyte protease,Fragmentin-2,Granzyme-2,Human lymphocyte protein,HLP,SECT,T-cell serine protease 1-3E,CGL1,CSPB,CTLA1,GRB
Accession No.	Uniprot ID:P10144
Uniprot	P10144
GeneID	3002;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:GZMB

Note: This product is for in vitro research use only