

Fibulin-1

Catalog No: #AP83525



Package Size: #AP83525-1 50ug #AP83525-2 100ug #AP83525-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Fibulin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	FIBL-1
Accession No.	Uniprot ID:O73775
Uniprot	O73775
GeneID	373979;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: FBLN1

Note: This product is for in vitro research use only