

Hepatocyte growth factor

Catalog No: #AP83527



Package Size: #AP83527-1 50ug #AP83527-2 100ug #AP83527-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Hepatocyte growth factor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Hepatopoietin-A, Scatter factor, SF
Accession No.	Uniprot ID:Q76BS1
Uniprot	Q76BS1
GeneID	282879;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: HGF

Note: This product is for in vitro research use only