

Heat-stable enterotoxin receptor

Catalog No: #AP83528



Package Size: #AP83528-1 50ug #AP83528-2 100ug #AP83528-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Heat-stable enterotoxin receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	STA receptor,Guanylyl cyclase C,GC-C,Intestinal guanylate cyclase,GUC2C
Accession No.	Uniprot ID:P55204
Uniprot	P55204
GeneID	397193;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:GUCY2C

Note: This product is for in vitro research use only