

Alpha-2-HS-glycoprotein

Catalog No: #AP83530



Package Size: #AP83530-1 50ug #AP83530-2 100ug #AP83530-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Alpha-2-HS-glycoprotein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Asialofetuin, Fetuin-A, FETUA
Accession No.	Uniprot ID:P12763
Uniprot	P12763
GeneID	280988;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: AHSG

Note: This product is for in vitro research use only