

Glutaredoxin-3

Catalog No: #AP83532



Package Size: #AP83532-1 50ug #AP83532-2 100ug #AP83532-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Glutaredoxin-3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	PKC-interacting cousin of thioredoxin,PICOT,PKC-theta-interacting protein,PKCq-interacting protein,Thioredoxin-like protein 2,PICOT,TXNL2,HUSSY-22
Accession No.	Uniprot ID:O76003
Uniprot	O76003
GeneID	10539;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:GLRX3

Note: This product is for in vitro research use only