

Platelet glycoprotein Ib alpha chain

Catalog No: #AP83533



Package Size: #AP83533-1 50ug #AP83533-2 100ug #AP83533-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Platelet glycoprotein Ib alpha chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	GP-Ib alpha,GPIb-alpha,GPIbA,Glycoprotein Ibalpha,CD42b
Accession No.	Uniprot ID:Q28256
Uniprot	Q28256
GeneID	403638;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:GP1BA

Note: This product is for in vitro research use only