

# Mast cell surface glycoprotein Gp49A

Catalog No: #AP83534



Package Size: #AP83534-1 50ug #AP83534-2 100ug #AP83534-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mast cell surface glycoprotein Gp49A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Gp49
Accession No.	Uniprot ID:Q61450
Uniprot	Q61450
GeneID	14727;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Gp49a

Note: This product is for in vitro research use only