

Glypican-2

Catalog No: #AP83538



Package Size: #AP83538-1 50ug #AP83538-2 100ug #AP83538-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Glypican-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Cerebroglycan,HSPG M13
Accession No.	Uniprot ID:P51653
Uniprot	P51653
GeneID	171517;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Gpc2

Note: This product is for in vitro research use only