

Pancreatic secretory granule membrane major glycoprotein GP2

Catalog No: #AP83539

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP83539-1 50ug #AP83539-2 100ug #AP83539-3 1mg

Description

Product Name	Pancreatic secretory granule membrane major glycoprotein GP2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Pancreatic zymogen granule membrane protein GP-2
Accession No.	Uniprot ID:P25291
Uniprot	P25291
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: GP2

Note: This product is for in vitro research use only