

Ferritin heavy chain

Catalog No: #AP83547



Package Size: #AP83547-1 50ug #AP83547-2 100ug #AP83547-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ferritin heavy chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Ferritin H subunit, Cell proliferation-inducing gene 15 protein, FTH, FTHL6, OK/SW-cl.84, PIG15
Accession No.	Uniprot ID:P02794
Uniprot	P02794
GeneID	2495;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: FTH1

Note: This product is for in vitro research use only