

Follicle-stimulating hormone receptor

Catalog No: #AP83548



Package Size: #AP83548-1 50ug #AP83548-2 100ug #AP83548-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Follicle-stimulating hormone receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	FSH-R,Follitropin receptor,LGR1
Accession No.	Uniprot ID:P23945
Uniprot	P23945
GeneID	2492;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FSHR

Note: This product is for in vitro research use only