

Alpha-ketoglutarate-dependent dioxygenase FTO

Catalog No: #AP83549



Package Size: #AP83549-1 50ug #AP83549-2 100ug #AP83549-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Alpha-ketoglutarate-dependent dioxygenase FTO
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Fat mass and obesity-associated protein,KIAA1752
Accession No.	Uniprot ID:Q9C0B1
Uniprot	Q9C0B1
GeneID	79068;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FTO

Note: This product is for in vitro research use only