

EGF-containing fibulin-like extracellular matrix protein

1

Catalog No: #AP83551

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP83551-1 50ug #AP83551-2 100ug #AP83551-3 1mg

Description

Product Name	EGF-containing fibulin-like extracellular matrix protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Fibulin-3,FIBL-3,Fbln3
Accession No.	Uniprot ID:Q8BPB5
Uniprot	Q8BPB5
GeneID	216616;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Efemp1

Note: This product is for in vitro research use only