

Fibroleukin

Catalog No: #AP83555



Package Size: #AP83555-1 50ug #AP83555-2 100ug #AP83555-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Fibroleukin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Fibrinogen-like protein 2,pT49
Accession No.	Uniprot ID:Q14314
Uniprot	Q14314
GeneID	10875;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FGL2

Note: This product is for in vitro research use only