

# Guanine nucleotide-binding protein G(I)/G(S)/G(T) subunit beta-3

Catalog No: #AP83556

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP83556-1 50ug #AP83556-2 100ug #AP83556-3 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Guanine nucleotide-binding protein G(I)/G(S)/G(T) subunit beta-3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Transducin beta chain 3
Accession No.	Uniprot ID:P52287
Uniprot	P52287
GeneID	60449;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Gnb3

Note: This product is for in vitro research use only