

Hyaluronan synthase 1

Catalog No: #AP83558



Package Size: #AP83558-1 50ug #AP83558-2 100ug #AP83558-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Hyaluronan synthase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Hyaluronate synthase 1,Hyaluronic acid synthase 1,HA synthase 1,Has
Accession No.	Uniprot ID:Q61647
Uniprot	Q61647
GeneID	15116;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Has1

Note: This product is for in vitro research use only