

Glutaminase kidney isoform, mitochondrial

Catalog No: #AP83568



Package Size: #AP83568-1 50ug #AP83568-2 100ug #AP83568-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Glutaminase kidney isoform, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	GLS,Gls1,Kiaa0838
Accession No.	Uniprot ID:D3Z7P3
Uniprot	D3Z7P3
GeneID	14660;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Gls

Note: This product is for in vitro research use only