

# Guanine nucleotide-binding protein G(s) subunit alpha isoforms XLas



Catalog No: #AP83571

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP83571-1 50ug #AP83571-2 100ug #AP83571-3 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Guanine nucleotide-binding protein G(s) subunit alpha isoforms XLas
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Alternative gene product encoded by XL-exon,Gnas1
Accession No.	Uniprot ID:Q6R0H6
Uniprot	Q6R0H6
GeneID	14683;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Gnas

Note: This product is for in vitro research use only