

## Hemoglobin subunit alpha-A

Catalog No: #AP83575



Package Size: #AP83575-1 50ug #AP83575-2 100ug #AP83575-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Hemoglobin subunit alpha-A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	Alpha-A-globin,Hemoglobin alpha-A chain
Accession No.	Uniprot ID:P01994
Uniprot	P01994
GeneID	416652;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:HBAA

Note: This product is for in vitro research use only