

# Guanine nucleotide-binding protein G(s) subunit alpha isoforms short

Catalog No: #AP83579

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP83579-1 50ug #AP83579-2 100ug #AP83579-3 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Guanine nucleotide-binding protein G(s) subunit alpha isoforms short
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Alternative gene product encoded by XL-exon,Gnas1
Accession No.	Uniprot ID:Q9Z213
Uniprot	Q9Z213
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Gnas

Note: This product is for in vitro research use only