

G-protein coupled estrogen receptor 1

Catalog No: #AP83585



Package Size: #AP83585-1 50ug #AP83585-2 100ug #AP83585-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	G-protein coupled estrogen receptor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Chemoattractant receptor-like 2,Flow-induced endothelial G-protein coupled receptor 1,FEG-1,G protein-coupled estrogen receptor 1,G-protein coupled receptor 30,GPCR-Br,IL8-related receptor DRY12,Lymphocyte-derived G-protein coupled receptor,LYGPR,Membrane estrogen receptor,mER,CEPR,CMKRL2,DRY12,GPER,GPR30
Accession No.	Uniprot ID:Q99527
Uniprot	Q99527
GeneID	2852;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:GPER1

Note: This product is for in vitro research use only