

Histone H2B type 2-E

Catalog No: #AP83590



Package Size: #AP83590-1 50ug #AP83590-2 100ug #AP83590-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Histone H2B type 2-E
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	H2b 613
Accession No.	Uniprot ID:Q64524
Uniprot	Q64524
GeneID	319190;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Hist2h2be

Note: This product is for in vitro research use only