

High mobility group protein HMGI-C

Catalog No: #AP83598



Package Size: #AP83598-1 50ug #AP83598-2 100ug #AP83598-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	High mobility group protein HMGI-C
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	High mobility group AT-hook protein 2,Hmgic
Accession No.	Uniprot ID:P52927
Uniprot	P52927
GeneID	15364;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Hmga2

Note: This product is for in vitro research use only