

Glucocorticoid receptor

Catalog No: #AP83602



Package Size: #AP83602-1 50ug #AP83602-2 100ug #AP83602-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Glucocorticoid receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	GR,Nuclear receptor subfamily 3 group C member 1,GRL
Accession No.	Uniprot ID:Q9N1U3
Uniprot	Q9N1U3
GeneID	396740;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:NR3C1

Note: This product is for in vitro research use only