

# Interleukin-4

Catalog No: #AP83607



Package Size: #AP83607-1 50ug #AP83607-2 100ug #AP83607-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Interleukin-4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	IL-4,B-cell IGG differentiation factor,B-cell growth factor 1,B-cell stimulatory factor 1,BSF-1,Lymphocyte stimulatory factor 1,II-4
Accession No.	Uniprot ID:P20096
Uniprot	P20096
GeneID	287287;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:IL4

Note: This product is for in vitro research use only