

Retinoic Acid Receptor α/β (Phospho-Ser96) Antibody

Catalog No: #12545



Package Size: #12545-1 50ul #12545-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Retinoic Acid Receptor α/β (Phospho-Ser96) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	Retinoic Acid Receptor α/β (Phospho-Ser96) Antibody detects endogenous levels of Retinoic Acid Receptor α/β only when phosphorylated at Retinoic AcSer96
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human Retinoic Acid Receptor α/β (Phospho-Ser96)
Target Name	Retinoic Acid Receptor α/β
Modification	Phospho
Other Names	RARA, NR1B1, PML-RAR, RAR-alpha, RARalpha, Retinoic acid receptor alpha, Pml-rara, RAR, RAR alpha, Retinoic acid receptor, alpha
Accession No.	Swiss-Prot#: P10276/P10826NCBI Gene ID: 5914/5915
Uniprot	P10276
GeneID	5914;
Target Species	human
Calculated MW	50kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000

Images



Western Blot analysis of lysates of MCF7, using primary antibody at 1:1000 dilution.

Note: This product is for in vitro research use only