

## Hydroxymethylglutaryl-CoA synthase, cytoplasmic

Catalog No: #AP83610



Package Size: #AP83610-1 50ug #AP83610-2 100ug #AP83610-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Hydroxymethylglutaryl-CoA synthase, cytoplasmic
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	HMG-CoA synthase,3-hydroxy-3-methylglutaryl coenzyme A synthase,Hmgcs
Accession No.	Uniprot ID:P17425
Uniprot	P17425
GeneID	29637;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Hmgcs1

Note: This product is for in vitro research use only