

Procollagen galactosyltransferase 1

Catalog No: #AP83611



Package Size: #AP83611-1 50ug #AP83611-2 100ug #AP83611-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Procollagen galactosyltransferase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Collagen beta(1-O)galactosyltransferase 1,Glycosyltransferase 25 family member 1,Hydroxylysine galactosyltransferase 1,GLT25D1
Accession No.	Uniprot ID:A5PK45
Uniprot	A5PK45
GeneID	513167;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:COLGALT1

Note: This product is for in vitro research use only