

Growth/differentiation factor 3

Catalog No: #AP83614



Package Size: #AP83614-1 50ug #AP83614-2 100ug #AP83614-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Growth/differentiation factor 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	GDF-3, UNQ222/PRO248
Accession No.	Uniprot ID:Q9NR23
Uniprot	Q9NR23
GeneID	9573;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: GDF3

Note: This product is for in vitro research use only