

Gamma-glutamyltranspeptidase 1

Catalog No: #AP83615



Package Size: #AP83615-1 50ug #AP83615-2 100ug #AP83615-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Gamma-glutamyltranspeptidase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	GGT 1, Gamma-glutamyltransferase 1, Glutathione hydrolase 1, Leukotriene-C4 hydrolase, CD224, Ggt, Ggtp
Accession No.	Uniprot ID:Q60928
Uniprot	Q60928
GeneID	14598;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Ggt1

Note: This product is for in vitro research use only