

# Gastrin/cholecystokinin type B receptor

Catalog No: #AP83619



Package Size: #AP83619-1 50ug #AP83619-2 100ug #AP83619-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Gastrin/cholecystokinin type B receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CCK-B receptor,Cholecystokinin-2 receptor,CCK2-R,CCKRB
Accession No.	Uniprot ID:P32239
Uniprot	P32239
GeneID	887;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:CCKBR

Note: This product is for in vitro research use only