

Low affinity immunoglobulin epsilon Fc receptor

Catalog No: #AP83621



Package Size: #AP83621-1 50ug #AP83621-2 100ug #AP83621-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Low affinity immunoglobulin epsilon Fc receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	BLAST-2,C-type lectin domain family 4 member J,Fc-epsilon-R1I,Immunoglobulin E-binding factor,Lymphocyte IgE receptor,CD23,CD23A,CLEC4J,FCE2,IGEBF
Accession No.	Uniprot ID:P06734
Uniprot	P06734
GeneID	2208;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FCER2

Note: This product is for in vitro research use only