

High affinity immunoglobulin epsilon receptor subunit alpha



Catalog No: #AP83622

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP83622-1 50ug #AP83622-2 100ug #AP83622-3 1mg

Description

Product Name	High affinity immunoglobulin epsilon receptor subunit alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Fc-epsilon RI-alpha,FcERI,IgE Fc receptor subunit alpha,Fce1a
Accession No.	Uniprot ID:P12371
Uniprot	P12371
GeneID	25047;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Fcer1a

Note: This product is for in vitro research use only