

Heparin cofactor 2

Catalog No: #AP83626



Package Size: #AP83626-1 50ug #AP83626-2 100ug #AP83626-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Heparin cofactor 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Heparin cofactor II,HC-II,Protease inhibitor leuserpin-2,Serpin D1,Hcf2,Hcii
Accession No.	Uniprot ID:P49182
Uniprot	P49182
GeneID	15160;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Serpind1

Note: This product is for in vitro research use only