

High mobility group protein B3

Catalog No: #AP83628



Package Size: #AP83628-1 50ug #AP83628-2 100ug #AP83628-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	High mobility group protein B3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	High mobility group protein 2a,HMG-2a,High mobility group protein 4,HMG-4,Hmg2a,Hmg4
Accession No.	Uniprot ID:O54879
Uniprot	O54879
GeneID	15354;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Hmgb3

Note: This product is for in vitro research use only