

Gastrokine-2

Catalog No: #AP83632



Package Size: #AP83632-1 50ug #AP83632-2 100ug #AP83632-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Gastrokine-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Blottin,Down-regulated in gastric cancer,Trefoil factor interactions(z) 1,BLOT,GDDR,TFIZ1,UNQ465/PRO813
Accession No.	Uniprot ID:Q86XP6
Uniprot	Q86XP6
GeneID	200504;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:GKN2

Note: This product is for in vitro research use only